



## NATURAL CONVECTION OVEN BONC-104

# NATURAL CONVECTION OVEN BONC-104

## NATURAL CONVECTION OVEN



- 304 stainless steel, mirror polishing processing, easy to clean and maintain.
- Natural convection with low noise.
- PID controller with over temperature alarm and timing function ensures precise and reliable control.

Option:

- Independent over-temperature alarm system ensures experiments running safely.
- RS 485 connector can connect computer to save the data via software.

## SPECIFICATIONS

Model	BONC-104
Electrical Requirement	AC220V 50HZ
Temperature Range	RT+10 ~ 200°C
Display Resolution	0.1°C
Temperature Stability	±0.5°C
Ambient Temperature	+5 ~ 35°C
Power Consumption	2000W
Chamber Volume	140L
Interior Dimension (WxDxH)mm	520x500x575
Exterior Dimension (WxDxH)mm	660x720x880
Shelves	3(pcs)
Packing Dimension	77*87*117
CBM	0.78
Net Weight	82.00
Gross weight	123
Alt Name	Natural Convection Oven

## APPLICATIONS

LED Microprocessor Controller (with timing function)

Provided for desiccation, torrefaction, wax-melting and sterilization in mining industry, laboratories and scientific research institutes.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)