



NATURAL CONVECTION OVEN BONC-310

NATURAL CONVECTION OVEN BONC-310

NATURAL CONVECTION OVEN

With easy to clean interiors and natural convection circulation system for efficient heating of the chamber, this oven can be relied on for general laboratory applications.



Stainless steel inner chamber without scroll rolled
Silicone rubber door seal
Strong sealing
Nice design

SPECIFICATIONS

| Model | BONC-310 |
|---|--------------------|
| Internal dimension/capacity (HxWxD) mm / liter | 415x370x345 |
| External dimension (HxWxD) mm | 570x750x515 |
| Package dimension (HxWxD) mm | 690x830x600 |
| Voltage (V) | 220 |
| Power (W) | 1000 |
| Temp. Range (°C) | 50~200 |
| Vacuum degree (Pa) | <133 |
| Temp. Precision (°C) | ±0.5 |
| Floating Range (°C) | ±1% |
| Weight (KG) (Wooden shelves) | 75 |
| Price (USD) | 460 |
| Note: Digital display intelligent temperature controller, the timing range is 0~999 minutes | |
| Alt Name | Vacuum drying oven |

APPLICATIONS

Biochemical, Pharmaceutical, Healthcare, Agricultural, Environmental protection, disinfection, Sterilization, Drying, Roasting, Laboratory, Research, engineering, Industry.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com