



HOLE TYPE OPERATION LIGHT BOPL-203

HOLE TYPE OPERATION LIGHT BOPL-203

HOLE-TYPE OR LAMP

Hole type operation light is suitable for small operating rooms. Labstac lights uses multiple different sources and by focusing the knob, the beam is illuminated from different surgery locations which reduces the shadow caused by medical staff head, hands and other parts caused. The illumination adjusted according to actual needs.



Color rendering Index (Ra) $100 \geq Ra \geq 85$

SPECIFICATIONS

Model	BOPL-203
Illumination	80,000LUX
Color temperature	4500K±500K
Color rendering Index	$100 \geq Ra \geq 85$
Depth of light beam	>500 cm
The Value of total irradiance	<6 mW/m ² lx
Total irradiance	<600 W/m ²
Temperature rise on the head of the patient under surgical operation	≤2°C
Rated Voltage of Bulb	24V
Rated Power of bulb	125 W
Input power	150 VA
Power Supply	220V±10% 50Hz
Alt Name	Hole-type OR lamp

APPLICATIONS

Medical unit, Surgery, Operation theatre



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com