



## INTEGRATED OPERATION LIGHT BOPL-302

# INTEGRATED OPERATION LIGHT BOPL-302

## INTEGRATED OR LAMP

Shadowless operation lamp is suitable for the requirement of various operation conditions and it is ideal operating lighting equipment for modern operation rooms. Based on trichromatic theory, the light is closer to natural light, which is suitable for medical workers to distinguish the color of vascular tissue, nerve ending and so on. Adapts US blue cold light coating technology, the optical filter could effectively filter ultraviolet rays from the spectrum.



Imported LED cold light source is adopted as operation lighting. As a genuine cold light source  
Streamline design  
The lamp can be flexible regulated 180 degrees  
Battery working (optional)

## SPECIFICATIONS

Model	BOPL-302
Illumination	80,000-160,000LX (at 1m distance)
Color temperature	4500K±500K
Temperature rise	≤12°C
Temperature rise on the head of the patient	≤2°C
Diameter of lamp holder	500 mm
Shift of major and stand by lamps	automatic
Depth of light beam	≥600mm
Facula diameter	φ100-φ280 mm
Brightness adjustment	8-stage continuous
Bulb Voltage	AC24V Bulb power
Input Power	200VA
Power Supply	AC220V±10% 50Hz
Alt Name	Integrated OR lamp

## APPLICATIONS

Urology, ENT, Gynecology examination, Surgery



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)