



## LED OPERATION LIGHT BOPL-404

# LED OPERATION LIGHT BOPL-404

## LED OR SHADOWLESS LAMP

Challenging tasks require special solutions. The innovative Labstac LED operation light give you benefits of cutting-edge LED technology in a uniquely effective manner. We focus on providing LED lights that are more color-accurate, more flexible, more economical. We offer innovative light and system solutions that can be retrofitted easily and without special preparation.



Stepless adjusting light source  
ONDAL arm, OSRAM bulb

## SPECIFICATIONS

Model	BOPL-404
Illumination	40000~160000 LUX/40000~160000 LUX Illumination at 1m on Endo mode- 10000 LUX
Color temperature	6700>Tc>3000K
Color rendering Index	100≥Ra≥85 Red color rendering Index R9≥97
Light field diameter	D10 180+40/180+40 mm
Depth of L1+L2 illumination	150 cm
Heat to light ratio	<6mW/m²lx
D50	110+35/110+35
Number of LED bulbs	84/84
Average bulb life (hours)	50000 hours
Ceiling height for installation	270-320 cm
Input power	150W/150W VA
Power Supply	100-240V
Alt Name	LED OR shadowless lamp

## APPLICATIONS

Dental implant surgery, Medical unit, Surgery, Operation theatre



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)