



## LED OPERATION LIGHT BOPL-408

# LED OPERATION LIGHT BOPL-408

## DOUBLE ARM LED OPERATING LAMP

Challenging tasks require special solutions. The innovative LED operation light gives you benefits of cutting-edge LED technology in a uniquely effective manner. We focus on providing LED lights that are more color-accurate, more flexible, more economical. We offer innovative light and system solutions that can be retrofitted easily and without special preparation.



## SPECIFICATIONS

| Model                                   | BOPL-408                      |
|---|-------------------------------|
| Illumination                            | 40000~120000 LuX              |
| Color Temperature (K)                   | 3000~5000                     |
| Color Rendering Index (Ra)              | ≥98                           |
| Special Color Rendering Index (R9)      | ≥98                           |
| Special Color Rendering Index (R13/R15) | ≥99                           |
| Diameter of Light Spot (mm)             | 120-350                       |
| Illumination Depth (mm)                 | ≥1500                         |
| Lamp Head Power (VA)                    | 180                           |
| LED Service Life (h)                    | ≥600000 Hour(s)               |
| Global Voltages                         | 100-240V 50/60 HZ             |
| Alt Name                                | Double Arm LED Operating Lamp |

## APPLICATIONS

Dental implant surgery, Medical unit, Surgery, Operation theatre.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)