



## SINGLE-HOLE LAMP BOPL-502

# SINGLE-HOLE LAMP BOPL-502

## SINGLE-HOLE LAMP

Single-hole lamp is suitable for urology, ENT, gynecology examination and general outpatient surgery partial auxiliary lighting. Easy to move with brakes and more suitable for cutting off power with a battery.



An imported LED cold light source is adopted  
Digital means are adopted for stepless regulation of LED brightness  
With electric focusing function  
The removable handle sheath can be disinfected at the high temperature of 135°C.  
Imported spring arm  
The imported switch power supply is adopted to control the voltage, making the work voltage permanently stable  
The no welding titanium alloy arm and the streamline-designed lamp cap better meet the design requirements of modern clean laminar flow operation rooms

## SPECIFICATIONS

Model	BOPL-502
Illumination	30000 LUX (at 1m distance)
Color temperature	4500K±500K
Color rendering Index	>85% Ra
Tilt left and right	180°
Tilt front and back	180°
Height adjustment distance	>400 mm
Input power	20+10% VA
Alt Name	Single-hole lamp

## APPLICATIONS

Dental implant surgery, Medical unit, Surgery, Operation theatre



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)