



MANUAL OPERATION TABLE BOPT-203

MANUAL OPERATION TABLE BOPT-203

MANUAL OR TABLE

Manual Operation table are suitable for various operation tables conducted by medical unit. The OR table and accessories are ergonomic and can be used across multiple disciplines. Various actions are controlled by hydraulic pneumatic system. 360° rotation helps every surgical requirement.



Various actions are controlled by hydraulic pneumatic system
All movements of the table are manipulated using the operating handle under the head plate
Imported Y type sealing ring is durable
The table surface adopts quality carbon plastic plate which is easily for conducting X-ray examination, equipped with coincidence rail
Both sides of the bed rail for the high quality 304 stainless steel
The base adopts imported carbon cloth which can be easily cleaned and sterilized
The mattress is wrapped with imported artificial leather. High-density memory sponge is formed in one piece, with no seam
Foot treadle brake device which is easily for moving operating table

SPECIFICATIONS

Model	BOPT-203
Table Length	2000±50 mm
Table Width	500±20 mm
Height Adjustment	(700~950) ±50 mm
Left tilt	≥15°
Right tilt	≥15°
Trendelenburg	≥20°
Reverse trendelenburg	≥15°
Head board upwards	≥30°
Head board downwards	≥90°
Backboard upwards	≥70°
Backboard downwards	≥15°
Waist board lifting	≥70°
Leg board downwards	≥90°
Leg board outwards	≥180°
Main Voltage AC	220V±22V, 50Hz±1Hz
Input Power VA	150±10%
Alt Name	Manual OR table

APPLICATIONS

Medical unit, Surgery, Operation theatre, cardiovascular, pediatric, gynecology, abdominal, gall bladder, plastic surgery and more



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com