



## MANUAL OPERATION TABLE BOPT-204

# MANUAL OPERATION TABLE BOPT-204

## MANUAL OR TABLE

Manual Operation table are suitable for various operation tables conducted by medical unit. The OR table and accessories are ergonomic and can be used across multiple disciplines. Various actions are controlled by hydraulic pneumatic system. 360° rotation helps every surgical requirement.



Imported LINAK motor various actions are controlled by remote controller  
The table, base seat and upright column covering are made of high quality 304 stainless steel, which is anti-corrosive and can be easily cleaned  
The mattress adopts high density memory sponge and is once formed, with no seam  
Whole stainless steel hidden auxiliary table with flexible extension and retraction  
Equipped with stainless steel sewage basin and other accessories  
Foot treadle brake device which is easy for moving obstetric table

## SPECIFICATIONS

Model	BOPT-204
Table Length	2000±50 mm
Table Width	600±20 mm
Height Adjustment	(700~950) ±50 mm
Working area length/width	520x500 mm
Trendelenburg/Rev-Trend	≥10°/≥12°
Backboard upwards	≥50°
Standard Accessory	A pair of arm support, a pair of leg support, a pair of foot pedal, a pair of arm rods, a stainless steel sewage basin, a handheld remote, controller
Main Voltage AC	220V±22V, 50Hz±1Hz
Input Power VA	150±10%
Alt Name	Manual OR table

## APPLICATIONS

Medical unit, Surgery, Operation theatre, cardiovascular, pediatric, gynecology, abdominal, gall bladder, plastic surgery and more



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)