



## ORTHOPEDICS TRACTOR RACK BOPT-403

# ORTHOPEDICS TRACTOR RACK BOPT-403

## ORTHOPEDICS TRACTOR RACK

Biolab's Orthopedics tractor rack is mainly applied to orthopedics surgical by using this tractor rack the patient's leg can be performed side traction, prone traction and open-up position of the body during the operation.



Adopt carbon fiber material, the perspective imaging is clear, light and strong  
The plug in design is convenient and quick to connect with table  
Equipped with a traction frame accessory car, which is convenient to move and save space

## SPECIFICATIONS

Model	BOPT-403
Table Length	1550±50 mm
Table Width	450±50 mm
Traction device horizontal rotary	0-180°
Traction distance	(0-200)+10%
Traction Device height adjustment	>190 mm
Alt Name	Orthopedics tractor rack

## APPLICATIONS

Surgery, Leg operation, Medical unit



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)