



VACUUM OVEN BOVA-7305

VACUUM OVEN BOVA-7305

VACUUM OVEN

Microprocessor controller (with timing function)

Vacuum oven is designed especially for drying of material which is thermo-sensitive or decomposes and oxidative easily. It can be filled with inert gases, which is especially for a rapid drying of some compound material.



- Microprocessor controller.
- Dual layer tempered glass door for clear observation.
- Minimum heating time 50% less than traditional vacuum oven.
- Interior chamber made from stainless steel with a mirror polished finish making for easy cleaning.
- Standard with vacuum pump.

SPECIFICATIONS

Model	BOVA-7305
Electrical Requirement	AC220V 50HZ
Controller	LED display
Power Consumption	1350W
Temperature Range	RT+10 ~ 200°C
Display Resolution	0.1
Temperature Stability	±1°C
Vacuum Degree	133Pa
Ambient Temperature	+5~40°C
Chamber Volume	90L
Interior Dimension (WxDxH) mm	450x450x450
Exterior Dimension (WxDxH) mm	610x590x1350
Shelves	2pcs (Independent temp. control) 4pcs
Chamber Material	304 Stainless steel
Packing Dimension	75*71*154
CBM	0.82
Net Weight	160.00
Gross weight	215
Alt Name	Vacuum Oven



MORE INFO

- Direct heating and control by shelves: Each shelf has a separate temperature sensor, can be independently displayed and controlled temperature in setting, monitoring and timing, more accurate and stable



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com