



## VACUUM OVEN BOVA-7402

# VACUUM OVEN BOVA-7402

## VACUUM OVEN

Programmability controller (with timing function)

Vacuum oven is designed especially for drying of material which is thermo-sensitive or decomposes and oxidative easily. It can be filled with inert gases, which is especially for a rapid drying of some compound material.



- Microprocessor controller.
- Dual layer tempered glass door for clear observation.
- Minimum heating time 50% less than traditional vacuum oven.
- Interior chamber made from stainless steel with a mirror polished finish making for easy cleaning.
- Standard with vacuum pump.
- Provides programmability for 7 periods and 9 steps for each period making for a total of 63 programmable step.

## SPECIFICATIONS

Model	BOVA-7402
Electrical Requirement	AC380V 50HZ
Controller	LCD Programmability controller
Power Consumption	3800W
Temperature Range	RT+10 ~ 200°C
Display Resolution	0.1
Temperature Stability	±1°C
Vacuum Degree	133Pa
Ambient Temperature	+5~40°C
Chamber Volume	431L
Interior Dimension (WxDxH) mm	630x810x845
Exterior Dimension (WxDxH) mm	1000x1040x1855
Shelves	4pcs (Independent temp. control)
Chamber Material	304 Stainless steel
Packing Dimension	114*120*203
CBM	2.78
Net Weight	508.00
Gross weight	656
Alt Name	Vacuum Oven





**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)