



## PIPETTE BULB BPIC-101

# PIPETTE BULB BPIC-101

RUBBER PIPETTE BULB



Our pinch and release valve system allows for one-hand operation  
Can accommodate pipettes up to 25mL  
Easy to use  
Made from natural rubber  
Chemical resistant rubber filler with 3 valves for suction  
Two valve system with a removable top valve for easy cleaning and one-handed use

## SPECIFICATIONS

Model	BPIC-101
Volume	2-25 ml
Color	Red
Unit	50
Alt Name	Rubber Pipette Bulb

## APPLICATIONS

Medical, Laboratory, Research, Chemistry, Biology, Liquid Handling Task, Institutes, Pharmaceutical, Industrial, Microbiology



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)