



PIPETTE PUMP BPIC-202

PIPETTE PUMP BPIC-202

PIPETTE PUMP

We are engaged in offering a wide range of Pipette Pumps that provide safe and accurate pipetting and dispensing of liquids. It can be rotated to dispense and draw up the liquid. Available in a variety of sizes and shapes.



- Economical
- Highly gas tightness
- Environment- friendly materials
- Lightweight and dexterous for operation
- Durable
- Recalibration Mechanism

SPECIFICATIONS

Model	BPIC-202
Volume	10ml
Colour	Green + White
Unit	1000
Alt Name	Pipette Pump

APPLICATIONS

Liquid Handling Task, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com