



PIPETTE STAND BPIC-302

PIPETTE STAND BPIC-302

PIPETTE STAND

Pipette stand is a holder manufactured to provide safe and convenient storage for pipettes. Modern design enables this stand to securely hold virtually every pipettes in the marketplace, making it truly 'universal'.



Linear type
Robust design
The freely rotating stand provides safe storage for all micropipettes
The stands are convenient to hold up 6 pipettes both of single and multi-channel.

SPECIFICATIONS

Model	BPIC-302
Description	For 6 pcs PS2 Series Fixed pipettes
Alt Name	Universal Pipette Stand

APPLICATIONS

Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task, Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com