



## PIPETTE TIPS BPIC-404

# PIPETTE TIPS BPIC-404

## PIPETTE TIPS

Pipette tips are essential when using a micropipette. These get attached to the shaft of the pipette and are usually ejected after use. It is an essential lab item for any test or experiment. They are available as racked or loose bulk tips, with or without a filter, sterile or non-sterile.



Made from medical grade Polypropylene  
Slippy surface with excellent design  
DNase-free  
Rnase-free  
Non-pyrogenic  
Filter and non-filter for choice  
Bag and rack are available to pack

## SPECIFICATIONS

Model	BPIC-404
Volume	1000 µl
Package	100 pcs/bag
Unit	10
Alt Name	Pipette Tips

## APPLICATIONS

Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)