



TIPS BOX BPIC-714

TIPS BOX BPIC-714

TIP BOX WITH ULTRA LOW RETENTION TIP

It offers quick and easy access to pipette tips, saving valuable time. Space-saving and environmentally compatible. Commonly used in liquid handling and molecular biology.



- Made from medical grade Polypropylene
- Slippy surface with excellent design
- DNase-free
- Rnase-free
- Non-pyrogenic
- Filter and non-filter for choice
- Bag and rack are available to pack

SPECIFICATIONS

Model	BPIC-714
Volume	300 µl
Unit	10
Alt Name	Tip box with ultra low retention tip

APPLICATIONS

Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com