



## BOTTLE TOP DISPENSER BPIP-101

# BOTTLE TOP DISPENSER BPIP-101

BOTTLE-TOP DISPENSER



Excellent chemical resistance  
Fully autoclavable at 121°C  
Adjustable control knob makes setting the volume easy  
Made of PTFE, FEP, BSG, PP  
Vapor pressure max. 500 mbar, viscosity max. 500 mm<sup>2</sup>/s, temperature max. 40°C,  
density max. 2.2 g/cm<sup>3</sup>.

## SPECIFICATIONS

Model	BPIP-101
Volume Range	2.5-25 ml
Subdivision	0.5 ml
Accuracy	±0.50%
Precision	0.10%
Volume	2.5-25 ml
Qty	1
Unit	pcs
Package	20 pcs/carton
Gross weight (kg)	14
Dimensions (m <sup>3</sup> )	0.595x0.43x0.275
Alt Name	Bottle-Top Dispenser



## ITEMS INCLUDED

No	Name
1	Adapt 45/40mm
2	Adapt 45/32mm
3	Adapt 45/38mm
4	Adapt 32/25mm
5	Adapt 32/28mm
6	Bottle (Brown, 1L)
7	Filling tube (m)
8	Discharge tube (m)



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)