



## LIQUID DISPENSER BPIP-204

# LIQUID DISPENSER BPIP-204

ADJUSTABLE GLASS-INJECTION DISPENSER



Availability for autoclave of the whole device  
A non-fragile, durable silicone tube made of special materials  
Convenience of use with high efficiency, and low repetitive error, capable of greatly reducing time and costs in research and experiment  
Solvent Resistance: general acid and alkali  
Material: PP

## SPECIFICATIONS

Model	BPIP-204
Bottle Volume	500 ml
Volume	0.5-5 ml
Test Volume	0.5 ml
Accuracy	±0.04
Qty	1
Unit	pcs
Package	10 pcs/carton
Gross weight (kg)	4.2
Dimensions (m <sup>3</sup> )	0.62x0.25x0.315
Alt Name	Adjustable Glass-Injection Dispenser amber glass



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)