



## PIPETTE CONTROLLER BPIP-411

# PIPETTE CONTROLLER BPIP-411

## MOTORIZED PIPETTE CONTROLLER

Pipette Controller offers an unmatched control of pipetting speeds from drop by drop or gravity dispense to fast liquid displacement. Designed for comfortable and uncomplicated pipetting even during long lasting sessions. It is light weight and easy to use in the laboratory and for research purposes.



- Chemical resistance
- Human engineering design
- Lightweight and Structure strong
- Lightweight and dexterous for operation
- A wide range volume for liquid controlling

## SPECIFICATIONS

|          |                             |
|----------|-----------------------------|
| Model    | BPIP-411                    |
| Volume   | 0.1-100 ml                  |
| Colour   | Dark blue                   |
| Alt Name | Motorize Pipette Controller |

## FEATURES

- Chemical resistance
- Human engineering design
- Lightweight and Structure strong
- Lightweight and dexterous for operation
- A wide range volume for liquid controlling



## APPLICATIONS

Liquid Handling Task, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)