



PIPETTE CONTROLLER BPIP-414

PIPETTE CONTROLLER BPIP-414

MOTORIZE PIPETTE CONTROLLER

Pipette Controller offers an unmatched control of pipetting speeds from drop by drop or gravity dispense to fast liquid displacement. Designed for comfortable and uncomplicated pipetting even during long lasting sessions. It is light weight and easy to use in the laboratory and for research purposes.



Chemical resistance
Human engineering design
Lightweight and Structure strong
Lightweight and dexterous for operation
A wide range volume for liquid controlling

SPECIFICATIONS

Model	BPIP-414
Volume	0.1-100 ml
Colour	Gray
Alt Name	Motorize Pipette Controller

FEATURES

Chemical resistance
Human engineering design
Lightweight and Structure strong
Lightweight and dexterous for operation
A wide range volume for liquid controlling



APPLICATIONS

Liquid Handling Task, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com