



SINGLE CHANNEL PIPETTE BPIP-501

SINGLE CHANNEL PIPETTE BPIP-501

FIXED VOLUME MICROPIPETTES

New design for smooth activation, excellent ergonomics, easy maintenance and fully autoclavable.



Precision calibration
Recalibration mechanism
Calibration adjustment with simple tools!

SPECIFICATIONS

Model	BPIP-501
Volume	1000 μ l
Accuracy $\leq$ %	0.8%
Accuracy $\pm$ μ l	8.0
CV $\leq$ %	0.20%
CV $\pm$ μ l	2.0
MOQ-20PCS	
Alt Name	Fixed Volume Single Channel Micropipette (Ejector 121°C autoclavable)

FEATURES

Precision calibration
Recalibration mechanism



Calibration adjustment with simple tools!





Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com