



SINGLE CHANNEL PIPETTE BPIP-502

SINGLE CHANNEL PIPETTE BPIP-502

STRAW MICROPIPETTE

The straw micropipette allows subsequent pipetting of 0.1 to 1.0 ml of the same liquid, this variable volume instrument provides greater cost and time efficiency. A universal metal nozzle accommodates straw tips (length 240 mm, diameter 4 mm), most adequate to aspirate from deep, narrow vials or stomacher bags, which is ideal for replacing glass and plastic pipettes for bacteriology applications.



Smooth activation, excellent ergonomics
Fully autoclavable 121°C / 250°F
Easy maintenance, cleaning and disinfection
Eject mechanism automatically
Digital volume showing

SPECIFICATIONS

Model	BPIP-502
Volume range	0.1-1 ml
Accuracy (%) - Min. vol.	< ± 2.0%
Accuracy (%) - Max. vol.	< ± 0.5%
Coefficient of variation (%) - Min. vol.	< 2.5%
Coefficient of variation (%) - Max. vol.	< 0.4%
MOQ-10PCS	
Alt Name	Straw Micropipette (121°C autoclavable)

FEATURES

Smooth activation, excellent ergonomics
Fully autoclavable 121°C / 250°F
Easy maintenance, cleaning and disinfection
Eject mechanism automatically
Digital volume showing



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com