



SINGLE CHANNEL PIPETTE BPIP-530

SINGLE CHANNEL PIPETTE BPIP-530

FIXED VOLUME MICROPIPETTES

New design for smooth activation, excellent ergonomics, easy maintenance and fully autoclavable.



Precision calibration
Recalibration mechanism
Calibration adjustment with simple tools!

SPECIFICATIONS

Model	BPIP-530
Volume	5 μ l
Accuracy $\leq$ %	2.0%
Accuracy $\pm$ μ l	0.1
CV $\leq$ %	0.80%
CV $\pm$ μ l	0.04
MOQ-20PCS	
Alt Name	Fixed Volume Single Channel Micropipette (121°C autoclavable)

FEATURES

Precision calibration
Recalibration mechanism



Calibration adjustment with simple tools!





Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com