



## SINGLE CHANNEL PIPETTE BPIP-535

# SINGLE CHANNEL PIPETTE BPIP-535

## FIXED VOLUME MICROPIPETTES

New design for smooth activation, excellent ergonomics, easy maintenance and fully autoclavable.



Precision calibration  
 Recalibration mechanism  
 Calibration adjustment with simple tools!

## SPECIFICATIONS

Model	BPIP-535
Volume	100 µl
Accuracy ≤ %	1.0%
Accuracy ± µl	1.6
CV ≤ %	0.25%
CV ± µl	0.25
MOQ-20PCS	
Alt Name	Fixed Volume Single Channel Micropipette (121°C autoclavable)

## FEATURES

Precision calibration  
 Recalibration mechanism



Calibration adjustment with simple tools!





**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)