



## SINGLE CHANNEL PIPETTE BPIP-538

# SINGLE CHANNEL PIPETTE BPIP-538

## FIXED VOLUME MICROPIPETTES

New design for smooth activation, excellent ergonomics, easy maintenance and fully autoclavable.



Precision calibration  
Recalibration mechanism  
Calibration adjustment with simple tools!

## SPECIFICATIONS

|                        |                                                               |
|------------------------|---------------------------------------------------------------|
| Model                  | BPIP-538                                                      |
| Volume                 | 500 $\mu$ l                                                   |
| Accuracy $\leq$ %      | 0.8%                                                          |
| Accuracy $\pm$ $\mu$ l | 8.0                                                           |
| CV $\leq$ %            | 0.20%                                                         |
| CV $\pm$ $\mu$ l       | 1.0                                                           |
| MOQ-20PCS              |                                                               |
| Alt Name               | Fixed Volume Single Channel Micropipette (121°C autoclavable) |

## FEATURES

Precision calibration  
Recalibration mechanism



Calibration adjustment with simple tools!





**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)