



SINGLE CHANNEL PIPETTE BPIP-543

SINGLE CHANNEL PIPETTE BPIP-543

FIXED VOLUME MICROPIPETTES

New design for smooth activation, excellent ergonomics, easy maintenance and fully autoclavable.



Precision calibration
Recalibration mechanism
Calibration adjustment with simple tools!

SPECIFICATIONS

Model	BPIP-543
Volume	25 μ l
Accuracy $\leq$ %	1.0%
Accuracy $\pm$ μ l	0.25
CV $\leq$ %	0.30%
CV $\pm$ μ l	0.075
MOQ-20PCS	
Alt Name	Fixed Volume Single Channel Micropipette (Ejector 121°C autoclavable)

FEATURES

Precision calibration
Recalibration mechanism



Calibration adjustment with simple tools!





Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com