



SINGLE CHANNEL PIPETTE BPIP-544

SINGLE CHANNEL PIPETTE BPIP-544

FIXED VOLUME MICROPIPETTES

New design for smooth activation, excellent ergonomics, easy maintenance and fully autoclavable.



Precision calibration
 Recalibration mechanism
 Calibration adjustment with simple tools!

SPECIFICATIONS

Model	BPIP-544
Volume	50 µl
Accuracy ≤ %	1.0%
Accuracy ± µl	0.5
CV ≤ %	0.30%
CV ± µl	0.15
MOQ-20PCS	
Alt Name	Fixed Volume Single Channel Micropipette (Ejector 121°C autoclavable)

FEATURES

Precision calibration
 Recalibration mechanism



Calibration adjustment with simple tools!





Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com