



SINGLE CHANNEL PIPETTE BPIP-545

SINGLE CHANNEL PIPETTE BPIP-545

FIXED VOLUME MICROPIPETTES

New design for smooth activation, excellent ergonomics, easy maintenance and fully autoclavable.



Precision calibration
Recalibration mechanism
Calibration adjustment with simple tools!

SPECIFICATIONS

Model	BPIP-545
Volume	100 μ l
Accuracy $\leq$ %	1.0%
Accuracy $\pm$ μ l	1.6
CV $\leq$ %	0.25%
CV $\pm$ μ l	0.25
MOQ-20PCS	
Alt Name	Fixed Volume Single Channel Micropipette (Ejector 121°C autoclavable)

FEATURES

Precision calibration
Recalibration mechanism



Calibration adjustment with simple tools!





Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com