



## SINGLE CHANNEL PIPETTE BPIP-548

# SINGLE CHANNEL PIPETTE BPIP-548

## FIXED VOLUME MICROPIPETTES

New design for smooth activation, excellent ergonomics, easy maintenance and fully autoclavable.



Precision calibration  
Recalibration mechanism  
Calibration adjustment with simple tools!

## SPECIFICATIONS

Model	BPIP-548
Volume	500 $\mu$ l
Accuracy $\leq$ %	0.8%
Accuracy $\pm$ $\mu$ l	8.0
CV $\leq$ %	0.20%
CV $\pm$ $\mu$ l	1.0
MOQ-20PCS	
Alt Name	Fixed Volume Single Channel Micropipette ( Ejector 121°C autoclavable)

## FEATURES

Precision calibration  
Recalibration mechanism



Calibration adjustment with simple tools!





**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)