



## PORTABLE REFRIGERATOR BPOR-100

# PORTABLE REFRIGERATOR BPOR-100

## PORTABLE REFRIGERATOR

The Portable Refrigerator is a versatile cooling solution designed for travel, camping, medical, and laboratory use, providing precise and reliable temperature control in a compact, durable design. Ideal for storing vaccines, food, and other temperature-sensitive items, this refrigerator combines portability with energy efficiency to ensure stable cooling in any environment.



## SPECIFICATIONS

Model	BPOR-100
Capacity	6L
Internal&External Materials	PP+PE
Insulation Material	PU
Temp. Resolution	0.1°C
Manufacturing Process	Blow molding+injection molding, burr-free surface
Cooling Source Material	Sodium polyacrylate
Cooling Retention Time	Maintains an internal temperature of 2~8°C for over 24 hours at an ambient temperature up to 43°C
Configuration	Refrigerator*1, phase-change ice box(400ml)*4, digital thermometer*1, shoulder strap*1
Internal Size(W*D*H)	240*145*180mm
External Size(W*D*H)	330*235*275mm
Package Size(W*D*H)	380*270*300mm
Net Weight	2.25kg
Gross Weight	2.7kg
Alt Name	Portable Refrigerator



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)