



PORTABLE REFRIGERATOR BPOR-108

PORTABLE REFRIGERATOR BPOR-108

PORTABLE REFRIGERATOR

The Portable Refrigerator is a versatile cooling solution designed for travel, camping, medical, and laboratory use, providing precise and reliable temperature control in a compact, durable design. Ideal for storing vaccines, food, and other temperature-sensitive items, this refrigerator combines portability with energy efficiency to ensure stable cooling in any environment.



SPECIFICATIONS

Model	BPOR-108
Capacity	85L
Internal&External Materials	PP+PE
Insulation Material	PU
Temp. Resolution	0.1°C
Manufacturing Process	Blow molding+injection molding, burr-free surface
Cooling Source Material	Sodium polyacrylate
Cooling Retention Time	Maintains an internal temperature of 2~8°C for over 24 hours at an ambient temperature up to 43°C
Configuration	Refrigerator*1, phase-change ice box(400ml)*15, digital thermometer*1
Internal Size(W*D*H)	580*465*320mm
External Size(W*D*H)	640*520*375mm
Package Size(W*D*H)	650*540*400mm
Net Weight	8kg
Gross Weight	9.3kg
Alt Name	Portable Refrigerator



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com