



TESTER PETROLEUM EQUIPMENT BPTL-241

TESTER PETROLEUM EQUIPMENT BPTL-241

WATER CONTENT TESTER

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.



The instrument is suitable to determine the water content of petroleum products and to determine the water content of lubricating grease. It adopts double-unit structure. It can do determination for two samples. The gripper is designed reasonably. The installation and dismantlement are convenient.

SPECIFICATIONS

Model	BPTL-241
Power supply	AC (220±10%) V, 50Hz
Total power consumption	≤ 2200W
Heating control	Continuously adjustable by a silicon knob, manual
Ambient temperature	(15-35) °C
Relative humidity	≤ 85%
Dimension	430 mm x 320 mm x 700 mm
Net weight	8.5 kg
Alt Name	Water content Tester

APPLICATIONS

Petroleum Industry, PVC Pipe Industry



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com