



## PARAFFIN TRIMMER BPTR-101

# PARAFFIN TRIMMER BPTR-101

## PARAFFIN TRIMMER



Full computer automated control.  
New heating element, fast heating, long life, energy saving.  
Temperature sensors are made in the USA, precise and reliable.  
Automatically saves the latest operational settings.  
LED displays the heating status.  
Permanent magnetic absorption L-shaped waste wax tank, built-in plastic film, easy to clean.

## SPECIFICATIONS

Model	BPTR-101
Temperature Range	Continuously adjustable within room temperature ~ 90°C
Working Area	230x160 mm
Working Voltage	AC 220V±10% 50Hz (standard model); AC110V±10% 60Hz
Power	150 W
Dimensions (LxWxH)	255x180x155 mm
Net Weight	3 kg
Alt Name	Paraffin Trimmer



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)