



PURE WATER SUPPLY SYSTEM BPWS-101

PURE WATER SUPPLY SYSTEM BPWS-101

PURE WATER SUPPLY SYSTEM

This system can be combined with any brand's water purification system to supply pure water for specific equipment and it is easy to install and maintain. It has Independent power control and it can be run automatically according to the liquid level of the tank.



Automatic microcomputer controlling system, LED real-time animation mode display.
 Running status is showed in the LED, such as flushing, producing water, full tank, water shortage, leakage and service.
 Power on self test, power reset, alarm when work more than 6 hours continuously, water shortage, leakage, low pressure and high pressure.
 3 procedure of the reverse osmosis membrane's self-flushing: power on, water shortage reset and work more than 2 hours continuously, extend the life of RO membrane.
 Bench top and floor stand(except for 45 series and built-in tank type), 2 kind installation method

SPECIFICATIONS

Model	BPWS-101
Shape & Volume of tank	Roundness: 100, 200 L Square:80,120, 200 L
Volume of Buffer Tank	12 L
Vertical self-priming height	2 meters
The Rated Stop Pressure	0.45 Mpa/0.4 Mpa
Max Pressure/Pump Lift	About 0.8 Mpa/80 meters
Alt Name	Pure Water Supply System

FEATURES

Automatic microcomputer controlling system, LED real-time animation mode display.
 Running status is showed in the LED, such as flushing, producing water, full tank, water shortage, leakage and service.
 Power on self test, power reset, alarm when work more than 6 hours continuously, water shortage, leakage, low pressure and high pressure.
 3 procedure of the reverse osmosis membrane's self-flushing: power on, water shortage reset and work more than 2 hours continuously, extend the life of RO membrane.
 Bench top and floor stand(except for 45 series and built-in tank type), 2 kind installation method
 High-strength shell with powder painting technics, achieve elegant appearance and meeting GLP standard
 Pretreatment cartridges, RO module, deionized cartridges, all designed to modularization independently. Easy to maintenance and replacement.
 Built-in 12 liters pressure tank (IT series), save lab space and easy to maintain.
 Different external tanks (optional) to meet every need and assure ample water-supply.
 Pipeline and fast-plug adaptor with NSF authorization, assure high quality ultrapure water.
 DOW's RO membrane, ensure stable operation and high desalinization rate.
 Precision polishing mixed resin cartridge, combine high pure water quality and low running cost.
 Portable TDS/conductivity test pen, testing feed water, RO water and deionized water's quality.

APPLICATIONS

Laboratory, Manufacturing, Reefkeeping, Aquarium, Laboratory, Research



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com