



HOT AIR STERILIZER BSHA-209

HOT AIR STERILIZER BSHA-209

PORTABLE PRESSURE STEAM STERILIZER

Durable, corrosion resistant and easy to clean stainless steel design in a wide range of capacities. Our hot air sterilizers are microprocessor controlled and perform reliably. A digital display with a broad observation window allows easy operation.



SPECIFICATIONS

Model	BSHA-209
Capacity	18 L
Function	Electric heated
Sterilization Temperature	126°C
Steam Exhausting Temperature	≤1°C
Working Pressure	0.14-0.16Mpa
Maximum Working Pressure	0.165 Mpa
Chamber Dimension	280x260 mm
Overall Dimension (LxWxH)	410x410x520 mm
Chamber Material	Stainless Steel
Weight	14/16 kg
Power	2 KW
Power Supply	220 V, 50 Hz
Alt Name	Portable Pressure Steam Sterilizer

APPLICATIONS

Medical, Research Institutions, Laboratory, Drinking Water Plant, Food industry, Hospital, Agriculture, Hospital, Scientific Research.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com