



## ORBITAL SHAKER BSOT-201

# ORBITAL SHAKER BSOT-201

## ORBITAL SHAKER

Orbital Shaker is designed with a compact appearance. It is used for any mixing applications in different industries. It is applied in different laboratories of molecular biology, microbiology, immunology, chemistry, etc.



DC brushless motor, safe, long life and maintenance-free. Operating is silent and stable.

Speed adjustable. Gentle shaking, applied in different industries.

Compact and clean designed, LED displays time and speed, streamlined fuselage, simple operation, easy clean.

Equipped with a timing function, which can be set arbitrarily within the range of 1min~ 100h. The time controller can automatically alarm, realizing unmanned operation, and can be set to run continuously from time to time.

Different platforms to load different conical flasks and other sample containers for optional choices.

## SPECIFICATIONS

Model	BSOT-201
Speed Range	50rpm~350rpm
Shaking Orbit	10mm
Time Range	1min~99h59min
Max. Loading Capacity	3kgs
Voltage	AC 100~230V, 50/60Hz
Power	50W
Fuse Protector	250V, 1A, $\Phi$ 5x20
Dimensions	W.284 x D.264 x H.108mm
Net weight	5.6kgs
Alt Name	Orbital Shaker

## ACCESSORIES FOR PURCHASE

No	Name
1	Adjustable Bar Platform
2	Flask Platform
3	Platform
4	Universal platform with springs
5	Flask Platform
6	Double Deck Platform
7	Sticky Pad



1



2



3



4



5



6



7



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)