



ORBITAL SHAKER BSOT-402

ORBITAL SHAKER BSOT-402

ORBITAL SHAKER

Orbital shakers are powerful variable speed shakers which provide efficient orbital motion. Speed and time are under microprocessor control. Shakers can be applicable in different labs: microbiology, chemistry, immunology, biochemistry and molecular biology.



Timer function, set the working time within 1 min-100 h.
Normal operation at 4~45 °C, safe and convenient.
A variety of pallets are available for different mixing operations.

SPECIFICATIONS

Model	BSOT-402
Shaking Speed	50-300 rpm
Orbit	10 mm
Max. capacity	3.0 kg
Timer range	1 min~99 h 59 min
Load capacity	See below table
Motor	DC brushless
Operating ambient	4~45 °C
Dimension (mm)	280x270x110 mm
Alt Name	Orbital Shaker

ACCESSORIES FOR PURCHASE

No	Name
1	flat platform with non-slip rubber mat
2	platform with clamps for flasks, 100-150 mL
3	platform with clamps for flasks, 250-300 mL
4	universal platform with adjustable bars for different types of flasks, bottles
5	platform with springs



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com