



PASS BOX BSPB-103

PASS BOX BSPB-103

PASS BOX



SPECIFICATIONS

Model	BSPB-103
Type	Pass Box
External Size(W*D*H)mm	860*770*780
Internal Size(W*D*H)mm	700*700*700
Material	304Stainless steel
Interlock	Electronical
UV Lamp	8W*1
HEPA Filter	/
Clean Air Speed	/
Power Supply	220V, 50Hz; 110V, 50/60Hz
Gross Weight (Kg)	90
Packing Size(W*D*H) mm	980*900*970
Alt Name	Pass Box

APPLICATIONS

Pass Box is used to transfer goods between rooms of different clean levels of clean room. It can reduce times of opening door, prevent the air convection between different rooms and minimize the pollution extent.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com