



ORBITAL SHAKER BSRP-104

ORBITAL SHAKER BSRP-104

ORBITAL SHAKER

Orbital shaker is a high speed shaker with direct brushless DC motor and microprocessor controlled technique. One touch knob for easy operation and can be changed to a different holder. It can be applicable in different laboratories, in microbiology, biotechnology, medical analysis, etc.



Maximum load capacity: 3 kg
Easy to operate with one touch knob
Variable speed between 50 to 300 rpm
DC brushless motor for gentle and powerful motion

SPECIFICATIONS

Model	BSRP-104
Speed Range	50-300 RPM
Timing Range	1 min ~ 99 h 59 min
Speed Increment	1 RPM
Loading	3 kg
Power	220 V / 120 V, 45 W
Dimensions (L x W x H) mm	280 x 270 x 110 mm
Weight	6.5 kg
MOQ-5PCS	
Alt Name	Orbital Shaker

ACCESSORIES FOR PURCHASE

No	Name
1	Flat Platform (MOQ-5PCS)
2	Platform with clamps for flasks (12 x 100 ml) (MOQ-5PCS)
3	Platform with clamps for flasks (6 x 250 ml) (MOQ-5PCS)
4	Universal platform with adjustable bars (MOQ-5PCS)



1



2



3



4



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com