



RECIPROCAL SHAKER BSRP-115

RECIPROCAL SHAKER BSRP-115

ROCKING SHAKER

The rocking shaker is a desktop oscillator that can be used directly in the greenhouse or in the incubator. It is widely used in the fixation of various electrophoresis gels, Coomassie blue staining, oscillating shaking during decolorization, silver nitrate staining Fix, dye, develop, etc.



The platform is made of stainless steel and natural rubber, which is acid & alkali resistant, non-deformation, easy-clean and spill-prevention.
Real-time display of speed and diminishing time.
The bearing is articulated, rotating, secure and no wear, suitable for long time use.
Compact designed, save space.
Simple and beautiful appearance, humanization design.
Adopt stepper motor, low noise.

SPECIFICATIONS

Model	BSRP-115
Speed Range	5~80rpm
Swing Cycloid	25mm (up and down) / 10°
Time Range	1min~99h59min
Max. Load	5.0kgs (double-deck platform)
Platform Dimension	L.280 x W.228mm
Voltage	AC 100~230V, 50/60 Hz
Power	50W
Fuse	250V, 1A, ϕ 5x20
Dimension	W.284 x D.264 x H.248mm
Net Weight	6.0kgs
Alt Name	Rocking Shaker

ITEMS INCLUDED

No	Name	Description
1	Double deck flat platform	with non-slip rubber mat



1



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com