



SINGLE BEAM UV VISIBLE SPECTROPHOTOMETER BSSBV-401

SINGLE BEAM UV VISIBLE SPECTROPHOTOMETER BSSBV-401

UV-VIS SPECTROPHOTOMETER



Standard scanning software can directly complete functions of Quantitative; Kinetics; Wavelength Scan; Multi Wavelength; DNA/Protein and Data processing.

Can establish calibration curves and implement associated tests. The instrument internal can be stored with 200 groups of data and 200 standard curves.

Suspended posture optical system design, strengthen and thicken the bottom plate to eliminate the vibration or transformation's impact on the optical system. Automatic wavelength calibration and automatic deviation repair.

Tungsten and Deuterium lamps can be changed easily, without adjustment.

Standard with PC software.

SPECIFICATIONS

Model	BSSBV-401
Wavelength Range	320-1100 nm
Bandwidth	2 nm
Wavelength Accuracy	±0.5 nm
Wavelength Reproducibility	≤0.2 nm
Photometric Accuracy	±0.3 %T
Photometric Repeatability	≤0.15 %T
Straylight	≤0.05 %T
Stability	±0.001A/h(500nm)
Baseline Flatness	±0.001 A
Noise	±0.0005 A
Photometric Range	0-200%T, -0.3-3A, 0-9999C
Wavelength setting mode	Automatic
Scanning speed	-
Output	USB Port
Printer port	Parallel Port
Display	LCD (128x64)
Light Source	Tungsten Lamp
Detector	Silicon Photodiode
Power	220V AC ±10%/50Hz or 110V AC / 60Hz
Dimension	460x380x180 mm
Weight	15 kg
Automatically bandwidth can be adjusted	-
Alt Name	Visible Spectrophotometer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com