



## BENCHTOP SHAKING INCUBATOR BSTH-109

# BENCHTOP SHAKING INCUBATOR BSTH-109

## MICROPLATE SHAKER

Shaker for microplates is a type of high speed shaker, which is with brushless DC motor and microprocessor controlled technique, specially for Elisa plates(96/384 wells, tissue culture plates (24/48/96 wells). The shaker is easy to operate with one touch knob and it can load 4 microplates.



- Continuous or timed operation
- Easy to operate with one touch knob
- Variable speed between 200 to 1,500 rpm
- DC brushless motor for gentle and powerful motion
- Orbital of 3 mm for efficient mixing in small sample vessels

## SPECIFICATIONS

Model	BSTH-109
Speed Range	200-1500 RPM
Timing Range	1 min ~ 99 h 59 min
Speed Increment	10 RPM
Mixing Orbit	3 mm
Capacity	4 microplate
Power	220 V / 120 V, 2 A
Dimensions (L x W x H) mm	280 x 270 x 110 mm
Weight	6.5 kg
MOQ-5PCS	
Alt Name	Microplate Shaker



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)