



BIOSAFETY TRANSPORT BOX BTBT-101

BIOSAFETY TRANSPORT BOX BTBT-101

BIOSAFETY TRANSPORT BOX



1. Real time temperature value display
2. Refrigeration and temperature control to ensure the consistency of samples
3. Easy to operate and could be used repeatedly
4. Meet the basic requirements of cold chain transportation and storage of vaccines

SPECIFICATIONS

Model	BTBT-101
Capacity	6L
Internal & External Materials	PP+PE
Insulation Material	PU
Temp. Resolution	0.1°C
Manufacturing Process	Blow molding + injection molding, burr-free surface
Cooling Source Material	Sodium polyacrylate
Cooling Retention Time	Maintains an internal temperature of 2~8°C for over 24 hours at an ambient temperature up to 43°C
Sponge Test Tube Holder	Standard: 20 holes, compatible with 10ml test tubes; Optional: 11 holes, compatible with 20ml test tubes
Configuration	Box*1, phase-change ice box(400ml)*4, absorbent cotton*1, shoulder strap*1, digital thermometer*1, sponge test tube holder*1
Internal Size(W*D*H)	240*145*180mm
External Size(W*D*H)	330*235*275mm
Package Size(W*D*H)	380*270*300mm
Net Weight	2.5kg
Gross Weight	3.2kg
Alt Name	Biosafety Transport Box

APPLICATIONS

It can be used for the transportation of biological products requiring physical insulation, such as un2814, un2900, un3373 biological samples, pathogenic microorganism(virus) species, blood, etc.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com