



TUBE ROLLER / ROTATOR BTUR-202

TUBE ROLLER / ROTATOR BTUR-202

LCD DIGITAL TUBE ROTATOR



LCD display of speed and time
Circular rotating action, gives gentle but effective mixing
Adjustable speed range from 10 to 70 rpm
Adjustable angle from 0 to 90°
Brushless motor ideal for extended use, and longer life
Wide variety of accessories are available

SPECIFICATIONS

Model	BTUR-202
Motor type	Brushless DC motor
Tilt angle	0-90°
Speed range	10-70rpm
Speed display	LCD
Timer display	LCD
Timer	1-1199min
Operation mode	Continuous / timer
Voltage, Frequency	100-240V, 50/60Hz
Power	40W
Dimension [DxWxH]	280 x 210 x 300mm
Permissible ambient temperature and humidity	5-40°C, 80%RH
Protection class	IP21
Alt Name	LCD Digital Tube Rotator



ACCESSORIES FOR PURCHASE

No	Name	Description
1	Disk Accessories Type-1	for 1.5mLx60 centrifuge tubes holder
2	Disk Accessories Type-2	for 15mLx16 centrifuge tubes holder
3	Disk Accessories Type-3	for 50mLx8 centrifuge tubes holder



1



2



3

APPLICATIONS

Laboratory Tube Rotators are ideal for Biochemistry, Molecular Biology and Clinical applications. Provides a gentle but effective mixing of biological samples in 1.5mL to 50mL micro tubes, used in a variety of applications including immune precipitations, prevention of blood coagulation, latex diagnostics etc.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com