



## TUBE ROLLER / ROTATOR BTUR-203

# TUBE ROLLER / ROTATOR BTUR-203

## ECONOMICAL ROTATOR



Provides gentle but effective rolling  
Adjustable speed range from 0 to 80 rpm

## SPECIFICATIONS

Model	BTUR-203
Motor type	DC motor
Speed range	0-80rpm
Speed display	Scale
Operation mode	Continuous
Voltage, Frequency	100-240V, 50/60Hz
Power	20W
Dimension [DxWxH]	150 x 530 x 190mm
Permissible ambient temperature and humidity	5-40°C, 80%RH
Protection class	IP21
Alt Name	Economical Rotator

## ITEMS INCLUDED

No	Name	Description
1	Plastic Clips	Disk/Rotisserie accessory, for 50mL, 15mL and 1.5mL



1

## APPLICATIONS

Horizontal Tube Rotator provides slow to vigorous mixing of your biological samples and solutions with advanced features you expect. Used in a variety of applications including immune precipitations, prevention of blood coagulation, latex diagnostic tests etc. Not recommended for extended use in incubators or cold rooms.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)