

TUBE ROLLER / ROTATOR BTUR-204

TUBE ROLLER / ROTATOR BTUR-204

HORIZONTAL LCD DIGITAL TUBE ROTATOR



LCD displays speed and time
Adjustable speed range of 10 to 70 rpm
Provides gentle but effective rolling and mixing
Can be used in cold rooms or incubators
Wide variety of Disk accessories and Rotisserie accessories available


SPECIFICATIONS

Model	BTUR-204
Motor type	Brushless DC motor
Speed range	10-70rpm
Speed display	LCD
Timer display	LCD
Timer	1-1199min
Operation mode	Continuous / timer
Voltage, Frequency	100-240V, 50/60Hz
Power	40W
Dimension [DxWxH]	220 x 510 x 260mm
Permissible ambient temperature and humidity	5-40°C, 80%RH
Protection class	IP21
Alt Name	Horizontal LCD Digital Tube Rotator



ACCESSORIES FOR PURCHASE

No	Name	Description
1	Disk support rods	4pcs in a group use with disk accessories
2	Rotisserie Accessory Type-1	for 1.5mLx48 centrifuge tubes held horizontally
3	Rotisserie Accessory Type-2	for 15mLx24 centrifuge tubes held horizontally
4	Rotisserie Accessory Type-3	for 50mLx24 centrifuge tubes held horizontally
5	Rotisserie Accessory Type-4	for 1.5mLx32 centrifuge tubes held vertically
6	Rotisserie Accessory Type-5	for 15mLx16 centrifuge tubes held vertically
7	Rotisserie Accessory Type-6	for 50mLx16 centrifuge tubes held vertically



APPLICATIONS

Laboratory Tube Rotators are ideal for Biochemistry, Molecular Biology and Clinical applications. Provides a gentle but effective mixing of biological samples in 1.5mL to 50mL micro tubes, used in a variety of applications including immune precipitations, prevention of blood coagulation, latex diagnostics etc.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com