



VORTEX MIXER BTUR-301

VORTEX MIXER BTUR-301

3D ROTATING MIXER

3D Rotating Mixer is a uniquely designed instrument which can mix samples by horizon, vertical and any angled rotation with a variety of sample tubes. It is widely used in molecular hybridization, blood and other samples that need to be mixed.



Unique combination of adjustable rotating frame.
A variety of rotation modes, from mild type to somersault type rotation mode.
Widely used in molecular biology, histo-chemistry, bio-chemistry, clinical applications, etc.
Provide three types of rotating frame.
It is positive and negative rotation. Just gently press the rotating frame when the instrument is working.
Oil-free bearings are built in both ends of the bracket.

SPECIFICATIONS

Model	BTUR-301
Speed Range	18rpm (Fixed Speed)
Reciprocating Swing Angle Range	/
Time Range	/
Voltage	AC 220V/AC 110V, 50/60Hz
Power	4W
Fuse	/
Dimension (WxDxH)	W.240 x D.145 x H.210mm
Net Weight	1kgs
Alt Name	3D Rotating Mixer

ITEMS INCLUDED

No	Name	Description
1	RH-A	5 centrifuge tubes of 10ml/15ml(Φ 15mm- Φ 16mm) +6 centrifuge tubes of 5ml /7 ml(Φ 12mm- Φ 13mm)
2	RH-B	3 centrifuge tubes of 50ml +18 centrifuge tubes of 0.5ml
3	RH-C	18 centrifuge tubes of 1.5ml / 2.0ml



1



2



3



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com