



TUBE ROLLER / ROTATOR BTUR-304

TUBE ROLLER / ROTATOR BTUR-304

3D GYRATORY ROCKER

3D Gyrotory Rocker is an economical 3D gyrotory rocker with perfect 3-dimensional motion for gentle but thorough mixing. It is an ideal choice for mixing small sample. It suits for blood collection tubes and centrifuge tubes, applied for all blot related applications such as gel staining and destaining. The concave pad can firmly keep the tubes.



Gentle but thorough mixing.
Suits for blood collection tubes and centrifuge tubes.
Fixed speed and tilt angle.
Safe for use in cold rooms and incubators.
Note: This model is available for 1.5mlx60, 3mlx30, 7mlx21, 15ml x17.
Also suitable for use in cold closet.

SPECIFICATIONS

Model	BTUR-304
Speed Range	20rpm (fixed)
Shaking Orbit	3D / 20 (fixed)
Max. Load	0.8kgs
Voltage	AC 220V/AC 110V, 50/60Hz
Power	45W
Dimension	W.240 x D.170 x H.150mm
Net Weight	1.6kgs
Alt Name	3D Gyrotory Rocker



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com