



## VORTEX MIXER BTUR-402

# VORTEX MIXER BTUR-402

## VERTICAL ROTATING MIXER

Vertical Rotating Mixer provides efficient but soft mixing to keep the sample in a suspension status such as blood mixing. It is suitable for preventing blood clots, latex diagnosis, immunoprecipitation and similar applications.



Smart appearance, solid structure, simple operation. LED displays time and speed. DC brushless motor, stable performance, low noise and long using life. 360°vertical circulation rotation, operating smoothly and steadily. Different tube holder frames for optional to suit for different kinds tubes. Oil-free bearings are built in both ends of the bracket.

## SPECIFICATIONS

Model	BTUR-402
Speed Range	10rpm~80rpm
Shaking Orbit	circle
Time Range	1min~99h59min
Voltage	AC 100~230V, 50/60Hz
Power	16W
Fuse	250V, 1A, $\Phi$ 5x20
Dimension	W.470 x D.190 x H.160mm
Net Weight	2.8kgs
Alt Name	Vertical Rotating Mixer



## ACCESSORIES FOR PURCHASE

No	Name	Description
1	Holder Frame	1.5ml Centrifuge x 64
2	Holder Frame	15ml Sharp Bottom Tube x 22
3	Holder Frame	50ml Sharp Bottom Tube x 16



1



2



3



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)